

A close-up photograph of an hourglass with a red liquid flowing through it. The top bulb is filled with the red liquid, and a stream is falling into the bottom bulb. The background is a light blue gradient. The hourglass is positioned in the upper half of the frame, with the bottom bulb partially obscured by a purple gradient that covers the lower half of the image.

STA[®]-Liatest[®] FDP


Diagnostic and monitoring of hypercoagulable states

- Quantitative determination of Fibrin and Fibrinogen Degradation Products
- Fully automated assay
- Fast turnaround time: 8 minutes
- Liquid and ready-to-use reagent

STA[®]-Liatest[®] FDP

Quantitative determination of Fibrin and Fibrinogen Degradation Products (FDP)

STA[®]-Liatest[®] FDP is a new quantitative and fully automated assay for FDP determination. This assay represents a significant step forward in diagnosis and monitoring of hypercoagulable states.




Clinical applications

- Disseminated Intravascular Coagulation (DIC)
- Sepsis and septic shock (bacteria, virus, ...)
- Obstetrical complications
- Intoxication (poisons, drugs, ...)
- Cancer
- Immunological disorders...

Biological study

A study was performed on 98 intensive care patients with septic shock and suspected DIC. All samples (D0, D1 and D3) were pooled.

Graph of comparison of FDP level between control and septic shock group population



Significant increase of FDP level in the septic shock population compared to the control group population ($P < 0.0001$ – Mann-Whitney test).

Main benefits

- **Quantitative results:** quantitative antigen assay by immuno-turbidimetric method
- **Specificity:** use of monoclonal antibodies specific for Fibrin and Fibrinogen Degradation Products
- **Fully automated assay**
- **Large linear working range:** 4 to 150 µg/mL
- From normal values (< 5 µg/mL) to high pathological levels
- **Practical:**
- Liquid and ready-to-use reagent
- Dedicated calibrator and controls

REAGENT	CALIBRATOR	CONTROLS
STA[®]-Liatest[®] FDP (Cat. Nr. 00649) • 6 x 5 mL buffer (R1) • 6 x 5 mL latex (R2)	STA[®]-FDP Calibrator (Cat. Nr. 00655) • 2 x 5 calibration levels x 1 mL	STA[®]-FDP Control (Cat. Nr. 00654) • 6 x 2 control levels x 1.5 mL
On-board stability: • 3 days on STA-R [®] and STA Compact [®] with STA [®] -Mini Reducer At 2-8°C after opening: • 12 days in its original capped vial	On-board stability: • 4 hours on STA-R [®] and STA Compact [®]	On-board stability: • 24 hours on STA-R [®] and STA Compact [®] with STA [®] -Micro Reducer

BIBLIOGRAPHY

1. Yanada M., Matsushita T., Suzuki M., Kiyoi H., Yamamoto K., Kinoshita T., Kojima T., Saito H., Naoe T. "Disseminated intravascular coagulation in acute leukemia: clinical and laboratory features at presentation". *European Journal of Haematology*, 2006, 77, 282-287
2. Levi M., Toh CH., Thachil J., Watson HG. "Guidelines for the diagnosis and management of disseminated intravascular coagulation". *British Journal of Haematology*, 2009, 145, N°1, 24-33

For further information, please contact:



Diagnostica Stago S.A.S
 RCS Nanterre B305 151 409
 9, rue des Frères Chausson
 92600 Asnières sur Seine (France)
 Ph. +33 (0)1 46 88 20 20
 Fax +33 (0)1 47 91 08 91
 webmaster@stago.com
 www.stago.com