






minicap sebia
flex piercing

Capillary
electrophoresis
based on a solid
experience




sebia


MINICAP Flex Piercing




PROTEINS



HbA1c



Hb



CDT

Substantial technology advancement


Cap piercing

Sample mixing prior to analysis


- Whole blood analysis
 - no pre-analytical preparation
 - sampling on capped primary tubes



Labor-saving, safety



- Optimal whole blood homogeneity by inversion at the sampling point
- No gradient of red blood cells.
- Accurate measurement



MINICAP Flex Piercing

- **Multi-parameter instrument / large test menu**

- Proteins, Immunotyping, HbA1c, Hemoglobin disorders and CDT

- **Analytical performances**

- Based on a solid experience of more than 8,000 capillary electrophoresis systems worldwide

- **Fully automated**

- **Cap piercing:** maximizes productivity
- **Positive ID:** full traceability from primary tube to final result

- **Throughput with high resolution**

- 20 tests/h (Proteins)
- 8 tests/h (HbA1c and Hb)
- 10 tests/h (CDT)

- **Specifications**

- Dimensions : 44 x 41.5 x 58 cm
- Weight : 32 kg




- **Main functions of the SQL Client / Server database**

- Networking of all Sebia instruments
 - within a laboratory
 - between distant sites
- Remote validation
- Unlimited data storage
- Web access, Extranet...

- **Interpretation software and results validation**

- Intuitive and user friendly interface
- Multicriteria search function
- Automatic recall of patient's history
- Triage and automatic validation of the results



sebia Benelux SCS / Comm.V.

Jan Olieslagerslaan 41
1800 Vilvoorde
Belgique / België
Tel. : +32 (0)2 702 64 64
Fax : +32 (0)2 702 64 60
e-mail : sebia.benelux@sebia.be

sebia Brasil Ltda

Edificio Baker Office Tower
Rua Barão do Triunfo, 73 Conjunto 74
Bairro Brooklin Paulista, São Paulo, CEP 04602-000
Brasil
Tel. : +55 11 3849 0148
Fax : +55 11 3841 9816
e-mail : sebia@sebia.com.br

sebia GmbH

Münsterfeldallee 6
36041 Fulda
Deutschland
Tel. : +49 (0)661 3 30 81
Fax : +49 (0)661 3 18 81
e-mail : sebia@sebia.de

sebia Hispania S.A.

C/ Sicilia n° 394
08025 Barcelona
España
Tel. : +34 93 208 15 52
Fax : +34 93 458 55 86
e-mail : sebia@sebia.es

sebia Inc.

1705 Corporate Drive, Suite 400
Norcross, Georgia 30093
U.S.A.
Tel. : +1 770 446 - 3707
Fax : +1 770 446 - 8511
e-mail : info@sebia-usa.com

sebia Italia S.r.l.

Via Antonio Meucci, 15/A
50012 Bagno a Ripoli (FI)
Italia
Tel. : +39 055 24851
Fax : +39 055 0982083
e-mail : info@sebia.it

sebia UK Ltd

River Court, Meadows Business Park
Station Approach, Blackwater,
Camberley, Surrey, GU17 9AB
United Kingdom
Tel. : +44 (0) 1276 600636
Fax : +44 (0) 1276 38827
e-mail : sales@sebia.co.uk

sebia Shanghai Rep. Office

Cross Tower, Room 2306-07
318 Fuzhou Road
Shanghai 200001
China
Tel. : +86 (21) 6350 1655
Fax : +86 (21) 6361 2011
e-mail : sebia@sebia.cn

sebia

Parc Technologique Léonard de Vinci
CP 8010 Lisses - 91008 EVRY Cedex
France
Tél. : +33 (0)1 69 89 80 80
Fax : +33 (0)1 69 89 78 78
e-mail : sebia@sebia.com

sebia